



DICOM CONFORMANCE STATEMENT FOR THE LODOX STATSCAN SYSTEM

Document no ASP-100-702

Revision no 1.1

Copyright © 2004 by Lodox Systems

All rights reserved. No part of this document may be reproduced, translated, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without prior written permission of the owner.

Note:

If this copy is no longer in use, return to sender.

APPROVAL PAGE

| NAME | SIGNATURE | DATE |
|-----------------------------------|-----------|-----------------|
| Carlos Sousa Manager: Software | | 20 October 2003 |

AUTHOR/S

| NAME | SIGNATURE | DATE |
|---|-----------|-----------------|
| Ben Wright Senior Software Engineer | | 20 October 2003 |

DOCUMENT INFORMATION

DISTRIBUTION LIST

| COPIES: | OWNER | COPIES: | OWNER |
|---------|-------------------------|--------------|-------------------------|
| MASTER | Product Data Management | USERS (2) | Product Data Management |
| 1 | | 11 | |
| 2 | | 12 | |
| 3 | | 13 | |
| 4 | | 14 | |
| 5 | | 15 | |
| 6 | | 16 | |
| 7 | | 17 | |
| 8 | | 18 | |
| 9 | | 19 | |
| 10 | | 20 | |

ELECTRONIC SOURCE DATA

| | |
|-----------------------|--------------------------------------|
| EFFECTIVE DATE | 20 October 2003 |
| FILE SIZE | 202 Kbytes |
| FILENAME | LODOX Statscan DICOM Conformance.doc |
| APPLICATION | MS Word 2002 |

CHANGE HISTORY

| Date | CR No | Revision | Description of change |
|-----------|-------|----------|-----------------------------|
| 20-Oct-03 | - | 1.0 | New document |
| 02-Feb-04 | | 1.1 | Merge of separate documents |
| | | | |
| | | | |
| | | | |
| | | | |

TABLE OF CONTENTS

| | | |
|-----------|---|-----------|
| 1. | SCOPE | 6 |
| 1.1 | Introduction | 6 |
| 1.2 | Purpose | 6 |
| 1.3 | Applicability | 6 |
| 2. | APPLICABLE DOCUMENTS AND REFERENCES..... | 7 |
| 2.1 | Applicable Documents | 7 |
| 2.2 | References | 7 |
| 3. | ACRONYMS, ABBREVIATIONS AND DEFINITIONS..... | 8 |
| 3.1 | Acronyms and Abbreviations | 8 |
| 4. | IMPLEMENTATION MODEL..... | 9 |
| 4.1 | Application Data Flow Diagrams | 9 |
| 5. | FUNCTIONAL DEFINITION OF AES..... | 10 |
| 5.1 | Modality Worklist SCU | 10 |
| 5.1.1 | Modality Worklist SCU AE Specification | 10 |
| 5.1.2 | Association Establishment | 10 |
| 5.1.3 | Association Initiation Policy..... | 11 |
| 5.1.4 | Association Acceptance Policy | 12 |
| 5.2 | Store SCU..... | 13 |
| 5.2.1 | Store SCU AE Specification..... | 13 |
| 5.2.2 | Association Establishment | 13 |
| 5.2.3 | Association Initiation Policy..... | 14 |
| 5.2.4 | Association Acceptance Policy | 14 |
| 5.3 | Print SCU | 15 |
| 5.3.1 | Print SCU AE Specification..... | 15 |
| 5.3.2 | Association Establishment | 15 |

| | | |
|-----------|---|-----------|
| 5.3.3 | Association Initiation Policy..... | 15 |
| 5.3.4 | Association Acceptance Policy | 18 |
| 6. | COMMUNICATION PROFILES..... | 19 |
| 6.1 | Supported Communication Stacks..... | 19 |
| 6.1.1 | OSI Stack..... | 19 |
| 6.1.2 | TCP/IP Stack | 19 |
| 7. | EXTENSIONS, SPECIALIZATIONS AND PRIVATIZATIONS | 20 |
| 7.1 | Standard Extensions..... | 20 |
| 7.1.1 | Standard Extension to DX Image for Storage SOP Class..... | 20 |

1. SCOPE

1.1 Introduction

This DICOM Conformance Statement specifies the behavior and functionality of software contained within the Lodox Statscan full-body digital X-ray scanner. It includes both the Statscan System operating software and the Diagnostic Viewing System (DVS) Software

1.2 Purpose

This DICOM Conformance Statement documents the conformance of the Lodox Statscan software with the Digital Imaging and Communications in Medicine (DICOM) standard. This document is essential in order to evaluate whether or not another DICOM compliant device can communicate with this software. This statement conforms to the recommended format as described in PS 3.2 of the DICOM standard.

1.3 Applicability

Lodox Statscan System Software.

2. APPLICABLE DOCUMENTS AND REFERENCES

2.1 Applicable Documents

2.2 References

DICOM (Digital Imaging and Communications in Medicine), NEMA, PS3.1 – 16, 2001.

3. ACRONYMS, ABBREVIATIONS AND DEFINITIONS

3.1 Acronyms and Abbreviations

| | |
|--------|--|
| ASCII | American Standard Code for Information Interchange |
| AE | Application Entity |
| ANSI | American National Standards Institute |
| DICOM | Digital Imaging and Communications in Medicine |
| DVS | Diagnostic Viewing System |
| DX | Digital X-ray |
| IE | Information Entity |
| IOD | Information Object Definition |
| ISO | International Standards Organization |
| MWL | Modality Worklist |
| NEMA | National Electrical Manufacturers Association |
| OSI | Open Systems Interconnection |
| PDU | Protocol Data Unit |
| SCP | Service Class Provider |
| SCU | Service Class User |
| SOP | Service Object Pair |
| TCP/IP | Transmission Control Protocol / Internet Protocol |
| UID | Unique Identifier |
| VM | Value Multiplicity |
| VR | Value Representation |

4. IMPLEMENTATION MODEL

The Statscan system is able to query a Modality Worklist (MWL) server in order to obtain patient demographic information from the MWL server. That is, the Statscan system acts as a DICOM Modality Worklist SCU.

The DVS software can send images to a DICOM-compatible storage system, as well as print images to a DICOM-compatible printer/laser camera. That is, the DVS software acts as a DICOM Storage SCU and a DICOM Print SCU. It also can act as a DICOM Verification SCU. A single application entity is shared for the storage and printing SCUs.

4.1 Application Data Flow Diagrams




Figure 1: Implementation Models

5. FUNCTIONAL DEFINITION OF AES

5.1 Modality Worklist SCU

The Modality Worklist SCU Application Entity is started when an operator uses the Worklist Query application to search for Patient Demographic information on a remote server.

The query request (C-Find) is transmitted over a single association. The results of the query (if any) are returned to the SCU using the same association. When the SCP returns a status of Complete, the association is released and the Modality Worklist SCU terminates. If the transmission fails for any reason, the association is aborted and the Modality Worklist SCU terminates.

5.1.1 Modality Worklist SCU AE Specification

This application entity provides conformance to the following DICOM SOP classes as SCU:

| SOP Class Name | SOP Class UID |
|--|------------------------|
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 |

5.1.2 Association Establishment

5.1.2.1 General

The PDU Length is configurable by the System Administrator. If no value is specified, the PDU Length defaults to 16384 bytes.

SOP Class extended negotiation is not supported.

5.1.2.2 Number of Associations

The Modality Worklist SCU will only propose a single Association at a time.

5.1.2.3 Asynchronous Nature

Asynchronous operation is not supported by this SCU.

5.1.2.4 Implementation Identifying Information

The following implementation information is valid:

| | |
|--------------------------------|-----------------------------|
| Implementation Version Name | OFFIS_DCMTK_352 |
| Implementation Class UID | 1.2.276.0.7230010.3.0.3.5.2 |
| DICOM Application Context Name | 1.2.840.10008.3.1.1.1 |

5.1.3 Association Initiation Policy

The application entity initiates an association with the selected remote Modality Worklist SCP. The calling application entity name is configurable by the system administrator. If the calling application entity name has not been configured it is set to LODOX. The called application entity name must be configured together with the presentation address to be used by using the Remote Application Entity configuration tool.

5.1.3.1 Associated Real-World Activity

The user starts the Worklist Query application, enters fields to search with, chooses a Remote AE to establish the Association with and clicks "Query".

5.1.3.2 Proposed Presentation Contexts

For each of the supported SOP classes, the software proposes as SCU a single presentation context containing the following transfer syntaxes, any of which is acceptable:

| Transfer Syntax | UID | Role | Extended Negotiation |
|---------------------------|---------------------|------|----------------------|
| Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU | None |

5.1.3.3 SOP Specific Conformance for Modality Worklist Management

In the Worklist Query application, the following fields are editable:

- Scheduled Station AE Title.
- Procedure Start Date – using the Date Picker controls.
- Modality – may be CR, DX or zero-length.
- Performing Physician – mapped to Scheduled Performing Physician's Name (0040,0006).
- Patient Last Name – mapped to Last Name Component of Patient Name (0010,0010).
- Patient First Name – mapped to First Name Component of Patient Name (0010,0010).
- MRN / Patient ID – mapped to Patient ID (0010,0020). If this field is present then the Patient Name value will be zero-length.
- Accession Number – mapped to Accession Number (0008,0050). If this field is present then the Patient Name value will be zero length.

The following attributes are set as described for this SOP class:

| Attribute Name | Tag | Match |
|---------------------------------------|-------------|--|
| Scheduled Station AE Title | (0040,0001) | Single value or zero-length |
| Scheduled Procedure Start Date | (0040,0002) | Single value, range match or zero-length |
| Scheduled Procedure Start Time | (0040,0003) | Zero-length |
| Modality | (0008,0060) | Single value or zero-length |
| Scheduled Performing Physician's Name | (0040,0006) | Single value or zero-length |

| | | |
|--------------------------------------|-------------|-----------------------------|
| Scheduled Procedure Step Description | (0040,0007) | Zero-length |
| Scheduled Station Name | (0040,0010) | Zero-length |
| Scheduled Procedure Step Location | (0040,0011) | Zero-length |
| Pre Medication | (0040,0012) | Zero-Length |
| Scheduled Procedure Step ID | (0040,0009) | Zero-length |
| Requested Contrast Agent | (0032,1070) | Zero-length |
| Scheduled Procedure Step Status | (0040,0020) | Zero-length |
| Comments on the Sch. Proc. Step | (0040,0400) | Zero-length |
| Requested Procedure ID | (0040,1001) | Zero-length |
| Requested Procedure Description | (0032,1060) | Zero-length |
| Study Instance UID | (0020,000D) | Zero-length |
| Requested Procedure Priority | (0040,1003) | Zero-length |
| Patient Transport Arrangements | (0040,1004) | Zero-length |
| Accession Number | (0008,0050) | Single value or zero-length |
| Requesting Physician | (0032,1032) | Zero-length |
| Placer Order Number | (0040,2016) | Zero-length |
| Filler Order Number | (0040,2017) | Zero-length |
| Imaging Service Request Reason | (0040,2001) | Zero-length |
| Patient Name | (0010,0010) | Single value or zero-length |
| Patient ID | (0010,0020) | Single value or zero length |
| Issuer of Patient ID | (0010,0021) | Zero-length |
| Patient's Sex | (0010,0040) | Zero-length |
| Patient's Birth Date | (0010,0030) | Zero-length |
| Patient's Birth Time | (0010,0032) | Zero-length |
| Patient's Weight | (0010,1030) | Zero-length |
| Patient's Size | (0010,1020) | Zero-length |
| Confidentiality Constraint | (0040,3001) | Zero-length |
| Medical Alerts | (0010,2000) | Zero-length |
| Contrast Allergies | (0010,2110) | Zero-length |
| Pregnancy Status | (0010,21C0) | Zero-length |
| Special Needs | (0038,0050) | Zero-length |
| Patient State | (0038,0500) | Zero-length |
| Admission ID | (0038,0010) | Zero-length |
| Issuer of Admission ID | (0038,0011) | Zero-length |
| Visit Status ID | (0038,0008) | Zero-length |
| Current Patient Location | (0038,0300) | Zero-length |
| Referring Physician's Name | (0008,0090) | Zero-length |

5.1.4 Association Acceptance Policy

Not applicable, as the Statscan cannot accept incoming Associations.

5.2 Store SCU

The Store SCU Application Entity is started when a user requests transmission of an image in the DVS database to a remote location.

A single image is transmitted over a single association. When the transmission of the image is complete, the association is released and the Store SCU terminates. If the transmission fails for any reason, the association is aborted and the Store SCU terminates.

5.2.1 Store SCU AE Specification

This application entity provides conformance to the following DICOM SOP classes as SCU:

| SOP Class Name | SOP Class UID |
|--|-----------------------------|
| Digital X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |

5.2.2 Association Establishment

5.2.2.1 General

The PDU Length is configurable by the System Administrator. If no value is specified, the PDU Length defaults to 16384 bytes.

SOP Class extended negotiation is not supported.

Although both the Digital X-ray Image Storage For Presentation and Computed Radiography Image Storage SOP classes are supported, the Store SCU will only use the Computed Radiography Image Storage SOP class if the SCP does not support the Digital X-ray Image Storage For Presentation SOP class. That is, the default behavior is for the DVS to act as SCU for the Digital X-ray Image Storage For Presentation SOP class. The Computed Radiography Image Storage SOP class should be considered as a “fallback” SOP class only.

5.2.2.2 Number of Associations

The Store SCU will only propose a single Association at a time.

5.2.2.3 Asynchronous Nature

Asynchronous operation is not supported by this SCU.

5.2.2.4 Implementation Identifying Information

The following implementation information is valid:

| | |
|-----------------------------|-----------------------------|
| Implementation Version Name | OFFIS_DCMTK_352 |
| Implementation Class UID | 1.2.276.0.7230010.3.0.3.5.2 |

5.2.3 Association Initiation Policy

The application entity initiates an association with the selected remote Storage SCP. The calling application entity name is set to LODOX. The called application entity name must be configured together with the presentation address to be used by using the Remote Application Entity configuration tool.

5.2.3.1 Associated Real-World Activity

The user selects the “Store” Option from the DICOM sub-menu, chooses a Remote AE to establish the Association with, and clicks “OK”.

5.2.3.2 Proposed Presentation Contexts

For each of the supported SOP classes, the software proposes as SCU a single presentation context containing the following transfer syntaxes, any of which is acceptable:

| Transfer Syntax | UID | Role | Extended Negotiation |
|---------------------------|---------------------|------|----------------------|
| Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU | None |

5.2.3.3 SOP Specific Conformance for DX and CR Image Storage Class

The following attributes of the Image Pixel Module are set as described:

| Attribute Name | Tag | Type | VR | VM | Comment |
|----------------------------|-------------|------|----|----|---------------------------------------|
| Samples Per Pixel | (0028,0002) | M/M | US | 1 | Value is 1 |
| Photometric Interpretation | (0028,0004) | M/M | CS | 1 | Value is MONOCHROME2 |
| Pixel Aspect Ratio | (0028,0034) | MC/M | | | Sent if pixel aspect ratio is not 1\1 |
| Bits Allocated | (0028,0100) | M/M | US | 1 | Value is 16 |
| Bits Stored | (0028,0101) | M/M | US | 1 | Value is 14 |
| High Bit | (0028,0102) | M/M | US | 1 | Value is 13 |
| Pixel Representation | (0028,0103) | M/M | US | 1 | Value is 0 |

5.2.4 Association Acceptance Policy

Not applicable, as the DVS software cannot accept incoming Associations. The DVS software as Storage SCP is not implemented.

5.3 Print SCU

The Print SCU Application Entity is started when a user requests printing of the current image to a remote DICOM printer.

The Print SCU prints a single Film Box as part of a print job over a single association. When the printing of the image is complete, the association is released and the Print SCU terminates. If the transmission to the printer fails for any reason, the association is aborted and the Print SCU terminates.

5.3.1 Print SCU AE Specification

This application entity provides conformance to the following DICOM SOP classes as SCU:

| SOP Class Name | SOP Class UID |
|---|-----------------------|
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 |

This application entity does not provide conformance to any SOP classes as SCP.

5.3.2 Association Establishment

5.3.2.1 General

The PDU Length is configurable by the System Administrator. If no value is specified, the PDU Length defaults to 16384 bytes.

SOP Class extended negotiation is not supported

5.3.2.2 Number of Associations

The Print SCU will only propose a single Association at a time.

5.3.2.3 Asynchronous Nature

Asynchronous operation is not supported by this SCU.

5.3.2.4 Implementation Identification Information

The following implementation information is valid:

| | |
|--------------------------------|-----------------------------|
| Implementation Version Name | OFFIS_DCMTK_352 |
| Implementation Class UID | 1.2.276.0.7230010.3.0.3.5.2 |
| DICOM Application Context Name | 1.2.840.10008.3.1.1.1 |

5.3.3 Association Initiation Policy

The application entity initiates an association with the selected remote Print SCP. The calling application entity name is set to LODOX. The called application entity name

must be configured together with the presentation address to be used by using the Remote Application Entity configuration tool.

5.3.3.1 Associated Real-World Activity

The user selects the “Print” Option from the DICOM sub-menu, chooses a Remote AE to establish the Association with, and clicks “OK”.

5.3.3.2 Proposed Presentation Contexts

For each of the supported SOP classes, the software proposes as SCU a single presentation context containing the following transfer syntaxes, any of which is acceptable:

| Transfer Syntax | UID | Role | Extended Negotiation |
|---------------------------|---------------------|-------------|-----------------------------|
| Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU | None |

5.3.3.3 SOP Specific Conformance for Printer SOP Class

Immediately after successful association negotiation, the Print SCU issues an N-GET-RQ message in order to retrieve the contents of the well-known Printer SOP Instance. The attribute identifier list element of the message remains empty, which means that the Print SCP is requested to transmit the contents of all attributes of the well-known Printer SOP Instance. If this request fails, the connection with the printer is released. Following retrieval of the well-known Printer SOP Instance, the association is released.

5.3.3.4 SOP Specific Conformance for Basic Film Session SOP Class

After retrieval of the well-known Printer SOP instance, a new association is established and, if this is successful, the Print SCU creates a Basic Film Session by issuing an N-CREATE-RQ message. The following elements are sent as part of the N-CREATE DIMSE Service Element:

| Attribute Name | Tag | Type | VR | VM | Comment |
|-----------------------|-------------|-------------|-----------|-----------|-----------------------|
| Number of Copies | (2000,0010) | U/M | IS | 1 | Default is 1 |
| Print Priority | (2000,0020) | U/M | CS | 1 | Default is HIGH |
| Medium Type | (2000,0030) | U/M | CS | 1 | Default is CLEAR FILM |
| Film Destination | (2000,0040) | U/M | CS | 1 | Default is MAGAZINE |
| Film Session Label | (2000,0050) | U/U | LO | 1 | As specified by user |

If creation of the Basic Film Session fails, Print SCU releases the association with the printer. After a successful completion of the print job, Print SCU uses an N-DELETE DIMSE Service Element to delete the Basic Film Session SOP instance before releasing the association. Other requests are never sent. In particular, Print SCU never sends an N-ACTION message on Basic Film Session level.

5.3.3.5 SOP Specific Conformance for Basic Film Box Class

After successful creation of the Basic Film Session SOP instance, the Print SCU creates a Basic Film Box. The following elements are sent as part of the N-CREATE DIMSE Service Element:

| Attribute Name | Tag | Type | VR | VM | Comment |
|----------------------------------|-------------|------|----|----|-------------------------|
| Image Display Format | (2010,0010) | M/M | ST | 1 | Fixed at 'STANDARD\1,1' |
| Referenced Film Session Sequence | (2010,0500) | M/M | SQ | 1 | |
| > Referenced SOP Class UID | (0008,1150) | M/M | UI | 1 | |
| > Referenced SOP Instance UID | (0008,1155) | M/M | UI | 1 | |
| Film Orientation | (2010,0040) | U/M | CS | 1 | Default is LANDSCAPE |
| Film Size ID | (2010,0050) | U/M | CS | 1 | Default is 14INX17IN |
| Magnification Type | (2010,0060) | U/M | CS | 1 | Default is REPLICATE |
| Trim | (2010,0140) | U/U | CS | 1 | Default is YES |

If creation of the Basic Film Box fails, the Print SCU releases the association with the printer. The Print SCU never creates more than a single Basic Film Box in the context of one association. The Print SCU uses an N-ACTION-RQ message to request processing of the print job from the Print SCP. After a successful completion of the print job, the Print SCU uses an N-DELETE-RQ message to delete the Basic Film Box SOP instance before deleting the Basic Film Session SOP instance and releasing the association. Other requests are never sent.

As part of creation of the Basic Film Box SOP Instance, the Print SCP creates an Instance of the Basic Grayscale Image Box SOP Class. The Print SCU uses an N-SET DIMSE Service Element to update the Basic Grayscale Image Box SOP Instance. The following elements are sent as part of the N-SET DIMSE Service Element:

| Attribute Name | Tag | Type | VR | VM | Comment |
|--------------------------------|-------------|------|----|----|---------------------------------------|
| Image Position | (2020,0010) | M/M | US | 1 | |
| Basic Grayscale Image Sequence | (2020,0110) | M/M | SQ | 1 | |
| >Samples Per Pixel | (0028,0002) | M/M | US | 1 | Value is 1 |
| >Photometric Interpretation | (0028,0004) | M/M | CS | 1 | Value is MONOCHROME2 |
| >Rows | (0028,0010) | M/M | US | 1 | |
| >Columns | (0028,0011) | M/M | US | 1 | |
| >Pixel Aspect Ratio | (0028,0034) | MC/M | | | Sent if pixel aspect ratio is not 1\1 |
| >Bits Allocated | (0028,0100) | M/M | US | 1 | Value is 8 |
| >Bits Stored | (0028,0101) | M/M | US | 1 | Value is 8 |
| >High Bit | (0028,0102) | M/M | US | 1 | Value is 7 |
| >Pixel Representation | (0028,0103) | M/M | US | 1 | Value is 0 |
| >Pixel Data | (7FE0,0010) | M/M | OW | 1 | |
| Polarity | (2020,0020) | U/M | CS | 1 | Default is NORMAL |

5.3.4 Association Acceptance Policy

Not applicable, as the DVS software cannot accept incoming Associations. The DVS software as Print SCP is not implemented.

6. COMMUNICATION PROFILES

6.1 Supported Communication Stacks

DICOM Upper Layer of TCP/IP is supported. The DICOM Toolkit used is the OFFIS DCMTK Toolkit, version 3.5.2

6.1.1 OSI Stack

Not supported.

6.1.2 TCP/IP Stack

The TCP/IP stack from the underlying Operating System is used. Port 104 is used by all SCUs.

7. EXTENSIONS, SPECIALIZATIONS AND PRIVATIZATIONS

7.1 Standard Extensions

7.1.1 Standard Extension to DX Image for Storage SOP Class

The following attributes from other SOP classes are added to the DX Image for Storage SOP Instance:

| Attribute Name | Tag | VR | VM | Comment |
|----------------|-------------|----|----|-------------------------|
| Scan Speed | (0018,1300) | DS | 1 | Taken from NM Image IOD |