

**Statscan®**  
**Versatility!**

 **lodox®**  
critical imaging technology

As simply stated by Frost & Sullivan:  
“One of most versatile DR systems in the market”


## RadMP System Upgrade

Lodox® has released the new **RadMP** system upgrade to extend Digital Radiographic (DR) examination functionality and multi-purpose flexibility on Statscan® Critical Imaging Systems.

In response to the radiographic exam needs of users for both trauma and standard Emergency Department patients and to further enhance the capabilities of Statscan as the most versatile DR system on the market, the new **RadMP** system upgrade is now available for new and existing Statscan® systems.

**RadMP** consists of a flexible 20" x 26" (50 x 67 cm) radiolucent imaging table that is fully adjustable and can be used in both horizontal and vertical planes. **RadMP** also includes an easy to adjust hydraulic patient imaging chair and pedestal stand for upright sitting and standing imaging purposes, ensuring accurate positioning, patient comfort and maximum safety.

The **RadMP** system offers allows the x-ray technologist to easily position seated or standing patients for upper-body extremities, erect chests/abdomens and also weight-bearing spine and joint procedures comfortably on sitting patients of all sizes and/or standing patients up to 5'2" (160cm) tall.



**Screen:** Light Weight 5.5lb (2.5kg), Image Area 20x26" (50x67cm), Washable, Skull angles to 45 Deg PA/AP, Removable Arm Support  
**Chair:** Total Weight 88lb (40kg), Max Load 440lb (200kg), Washable, Hydraulic lift 10-236" (25-600cm), Platform Lift 5-18" (13-45cm)  
Removable Seat 17x21" (45x55cm), Lockable castors, Foot Pedals, Lift /Lower/Lock

Emergency Care • Trauma Radiography • General Radiography

Low Dose • Fast Digital X-Ray

lodox systems na, llc, 108 east lake st. south lyon, mi 48178 tel: (248) 4468480/3 fax: (248) 4370769 u.s. toll free: 886-61-lodox (56369)  
lodox systems (pty) ltd. p.o. box 651807, benmore, 2010, south africa tel: +27 11 4449118 fax: +27 11 2620085 [www.lodox.com](http://www.lodox.com)